



A Brand Division of American Tooth Industries

1200 Stellar Drive • Oxnard, CA 93033-2404
(805) 487-9868 • Toll Free (800) 235-4639 • FAX: (805) 483-8482
E-mail: info@americantooth.com • Web Site: www.americantooth.com



Justi Miracleflask™

The Most Revolutionary Curing System for Full and Partial Denture Restorations and Relines.

For Use With Microwave Oven - Exceptional Quality in Denture Production:

A Complete System For Curing Dentures and Partials. Special materials and special techniques are not required. Removal of wax and curing of Justi Denture Base are easily performed. Microwave processed resin provides a superior fitting denture.

Justi Miracleflask™ and the Microwave Oven:

The Justi Miracleflask™ is the most revolutionary curing system for dentures, partials and relines. It is made of a fiber reinforced plastic material that is especially designed for use in a microwave oven. It may also be used with the conventional water bath heat curing method. It will withstand extreme pressure under an hydraulic or hand press up to 3,000 psi. Compared with the conventional method, the microwave radiation technique using the Justi Miracleflask™ gives a more complete curing of resin. Unpolymerized residual monomer is greatly reduced producing a final resin product with greater hardness. Cooling time and method are the same as in the conventional method. When Justi fibered denture base material is used with the Miracleflask™, our tests show superior results. There are many advantages to using the Justi Miracleflask™.

Major Attractions of the Miracleflask™.

- More accurate and dimensionally stable prosthesis.
Uniform curing of the acrylic resin.
Rapid process
Complete and total curing of the acrylic resin with no significant free monomer residual.
Clean removal method without damage to flask.
Reduced energy costs.
Expedient service.

Specifications:

- Material is fiber reinforced plastic (FRP).
Thermal deformation temperature is 482°F (250°C).
Compression strength is 18 kig/mm².
Weight is 900 g.
Dimensions are 14 cm (W) X 12 cm (L) X 7 cm (H).



Other reliable Justi Products:

- Harden Plastic Teeth
Justi Imperial®
Anteriors & Posteriors
Justi Blend®
Anteriors & Posteriors
Dymon-Hue® HPT
Anteriors & Posteriors
Justi® Imperial® Plastic Facings
Pyroplast Veneering Material
Justi Tray Material
Justi Milling and Polishing Cream
Justi Micro-Link Denture Base
Material-High Impact
Justi Repair Material
Justi Opaquer
Wax Solvent No. 20
Namilon Crown and Bridge
Material
Justi Cold-Cure Ortho Resin-
High Impact
Justi-Pakto Crown & Bridge Resin
and many more...

Packaging:

Table with 2 columns: Item # and Description. Includes items like 067-44-0001 (Complete Kit), 067-44-0002 (Intro Kit), and 067-45-0002 (Microwave Oven).

Replacement Parts

Table with 2 columns: Item # and Description. Lists replacement parts like Nuts & Bolts, Miracleflask Top/Middle/Bottom Sections, Pressure Plate, Tooth Wrench, and Wooden Hammer.

\*If shade is not indicated - Regular Fibered Denture Base will be included.

Products for Better Dentistry

Additional questions, to place your order, or for the Dealer nearest you, please call the Tooth Express Order Desk at 800-628-1437.