




"SIMPLE INVENTIONS LAST FOREVER"

# JUSTI MIRACLEFLASK™

The Most Revolutionary Curing System for Full and Partial Denture Restorations and Relines.



## For Use With Microwave Oven - Exceptional Quality in Denture Production:

### A Complete System For Curing Dentures and Partials.

Special materials and special techniques are not required. Removal of wax and curing of Justi Denture Base are easily performed. Microwave processed resin provides a superior fitting denture.

### Justi Miracleflask™ and the Microwave Oven:

The **Justi Miracleflask™** is the most revolutionary curing system for dentures, partials and relines. It is made of a fiber reinforced plastic material that is especially designed for use in a microwave oven. It may also be used with the conventional water bath heat curing method. It will withstand extreme pressure under a hydraulic or hand press up to 3,000 psi. Compared with the conventional method, the microwave radiation technique using the **Justi Miracleflask™** gives a more complete curing of resin. Unpolymerized residual monomer is greatly reduced producing a final resin product with greater hardness. Cooling time and method are the same as in the conventional method. When Justi fibered denture base material is used with the **Miracleflask™**, our tests show superior results. There are many advantages to using the **Justi Miracleflask™**.

### Major Attractions of the Justi Miracleflask™:

- More accurate and dimensionally stable prosthesis.
- Uniform curing of the acrylic resin.
- Rapid process.
- Complete and total curing of the acrylic resin. no significant free monomer residual.
- Clean removal method without damage to the flask.
- Reduced energy costs.
- Expedient service.

### Specifications:

- Material is fiber reinforced plastic (FRP).
- Thermal deformation temperature is 482°F (250°C)
- Compression strength is 18 kg/mm<sup>2</sup>
- Weight is 900 g.
- Dimensions are 14 cm (W) X 12 cm (L) X 7 cm (H).

### Packaging:

Item #	Description
067-44-0001	Justi Miracleflask Complete Kit Contents: 2 Miracleflasks, 6 Nuts & Bolts, 1 Wood Mallet, 1 Wrench, 1 lb. Denture Base Powder (*specify shade) & 8 oz. Liquid.
067-44-0002	Justi Miracleflask Intro Kit (Flask, Nut, Bolt & Wrench - 1 ea.)
067-44-0010	Replacement Nuts & Bolts (3 Sets)
067-44-0013	Miracleflask Top Section
067-44-0014	Miracleflask Middle Section
067-44-0015	Miracleflask Bottom Section
067-44-0016	Miracleflask Pressure Plate
067-44-0027	Tooth Wrench
067-44-0031	Justi Miracleflask Hammer
067-45-0002	Microwave Oven

\*If shade is not indicated - Regular Fibered Denture Base will be included.

### OTHER RELIABLE JUSTI PRODUCTS:

Hardened Plastic Teeth: Justi Imperial® Anterior & Posterior Justi Blend® Anterior Dymon-Hue® HPT Anterior & Posterior	Pyroplast Veneering Material Namilon Crown and Bridge Resin Justi-Pakto Crown and Bridge Resin Justi Rainbow Liquid Justi Repair Material Justi tray Material Wax Solvent No. 20 Milling and Polishing Cream Justi Opaquer Justi Denture Base Material Justi Cold-Cure Ortho Resin - High Impact and many more...
Porcelana Low Fusing, Synthetic Porcelain	
Justi Imperial Plastic Facings	

### READ & PASS ALONG TO

\_\_\_\_\_

\_\_\_\_\_

Then file for reference

### Distributed by:

### Developed, Designed and Manufactured by:

## Justi Products

### A Division of American Tooth Industries

1200 Stellar Drive • Oxnard, CA 93033-2404

(805) 487-9868 • (800) 235-4639 • Fax (805) 483-8482

E-mail: info@americantooth.com • Web site: www.americantooth.com