



# Justi Products

A Brand Division of American Tooth Industries

1200 Stellar Drive • Oxnard, CA 93033-2404

(805) 487-9868 • (800) 235-4639 • FAX: (805) 483-8482

E-mail: info@americantooth • Web site: www.americantooth.com

## Imperial® & Imperial Vivo® Block Teeth Order Form in Imperial, Bleached and Imperial Vivo Shades

Order # \_\_\_\_\_ Cust. # \_\_\_\_\_

Bill To: \_\_\_\_\_

Address \_\_\_\_\_

City, State, Zip \_\_\_\_\_

Phone \_\_\_\_\_ Fax \_\_\_\_\_

Date \_\_\_\_\_

Ship To \_\_\_\_\_

Address \_\_\_\_\_

City, State, Zip \_\_\_\_\_

Phone \_\_\_\_\_ Fax \_\_\_\_\_

Molds	10 Imperial Shades									
	20	25	45	59	62	65	66	67	69	77
<b>30M 084-50-0701U UPPER (1 x 8)</b>										
<b>30M 084-50-0701L LOWER (1 x 8)</b>										

4 Bleached Shades			
17	18	19	21

13 Imperial Vivo Shades												
A1	A2	A3	A3.5	A4	B1	B2	B3	C1	C2	C4	D2	D3

Grp/Item #

Upper: 084-50-0701US Total Uppers \_\_\_\_\_ X \$ \_\_\_\_\_ = \$ \_\_\_\_\_

Lower: 084-50-0701LS Total Loweres \_\_\_\_\_ X \$ \_\_\_\_\_ = \$ \_\_\_\_\_

S = Shade (e.g. 084-50-0701U59)

Total Sets \_\_\_\_\_ GRAND TOTAL \$ \_\_\_\_\_


Signed \_\_\_\_\_

Comments: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_



### Justi® Imperial® and Justi Imperial Vivo® Block Teeth (1 x 8)

Made from highly cross-linked acrylic material, the block teeth are 4 blended posterior teeth joined together on a flat plane. The zero degree is a time saving factor in articulation and mounting of the blocks in a denture wax up. The cusplless carvings eliminate much of the lateral stress during use.

Available in 10 Imperial shades, 4 Bleached shades & 13 Imperial Vivo shades.