



# Flexible Clear Disc



AMERICAN  
TOOTH  
INDUSTRIES

## For Invisible Partial Denture


The Flexible Disc is intended for aesthetically superior semi-flexible removable partial denture frameworks, which blend perfectly with the patient's gums, providing better function, aesthetics, comfort, and durability.

### Indications:

Transparent partial denture frameworks, night guards, clasps, and other dental appliances processed by CAD/CAM milling machines for lightweight, aesthetic non-metallic restorations.

### Benefits for Labs:

- 100% digital for accurate fit
- Eliminates RPD in the wax-up stage, investing, in casting and finishing.
- Fast, clean, precise, and cost-effective
- Easy finishing and polishing with minimal effort
- Mill with the same procedure with PMMA, both wet and dry mill



### Benefits for Patients:

- Lightweight, durable and precise
- Comfortable and easy adaptation in the mouth
- Biocompatible and non-allergenic
- Clasps are flexible and tooth-bearing, virtually unbreakable.
- Occlusal rest reduce tissue-borne pressure to reduce sore sport and offer excellent retention

Part Number: FLD9820-CL

Dimension: 20mm x98.5mm (2 steps)

Shade: Transparent

For any questions please contact us at (805)487-9868  
or visit our website at [www.americanntooth.com](http://www.americanntooth.com)