

Justi

Hollywood Teeth[®]

IMPERIAL JUSTI BLEND VIFORM

Mold Chart of
*Lingualized Occlusion
Combos[™]*



MED/CERT
ISO 13485:2013

CE
0482

NEW PRODUCT FOR 2008



Hollywood Teeth Line Lingualized Occlusion Combos™


The Professional's Choice

Together with experts on lingualized occlusion techniques we have put together a selection of our posterior brands which are suitable to be used in manufacturing dentures, overdentures, and selective partials. Since the Imperial line of posterior teeth is made of a large selection of sizes we were able to create 9 different molds for lingualized occlusion. For Justi Blend we present 6 combinations, and for Viform only 2 combinations are offered. All these brands are made of a very hard plastic with a very dense cross linking formula. Justi guarantees superior wear ability for all of our posterior brands.

IMPERIAL® Lingualized Occlusion Combos™

These combinations have been selected by experts to achieve the proper lingualized tilt for more efficient mastication. Our combinations feature the popular 33° upper vs. 10° lower molds and can be easily setup and mounted whether the application is for a denture, overdenture, or partial.

Imperial teeth are offered in our standard Eterna Age Shades system from #17 to #82. Many shades offered in the Eterna system are comparable to shades in other very popular shade guide systems. See the Imperial mold chart "Shade Comparison Chart" for guidance.



HOW TO READ THE TOOTH CARDS

A - Straight line measurement of 4 upper posteriors

B - Straight line measurement of 4 lower posteriors

Mold Number	30L	
Cusps Degrees	33°	
Total Width Upper (A)	29.9	Profile of 1 st Premolar
Mold Number	30X	
Cusps Degrees	10°	
Total Width Lower (B)	33.0	

Tooth measurements are within plus or minus 5% of stated value.

Justi Teeth Products



The Most Comprehensive Guarantee in the Industry

The Justi Company goes to great lengths to ensure the finest quality teeth. As a symbol of our commitment to quality, we provide the most comprehensive guarantee in the industry. Our teeth meet and exceed ANSI/ADA Specification 15, published in 1999 identical to ISO 3336, BS 3990, NF S91-213, DIN 13914, International Standards. For the certification of synthetic resin teeth to be used in dentures, partials and overdentures for implant technology. Justi teeth are stain resistant, wear and solvent resistant, color stable, fluorescent and photogenic. Justi teeth are warranted against excessive wear, breakage or discolorization for the life of your denture. This warranty does not extend to the denture base which supports the teeth nor does it cover the laboratory service fee or the dentist fee. Verification of laboratory records is required. A certificate of authenticity and warranty is included with each combination.

VIFORM® *Lingualized Occlusion Combos™*

The ViForm Lingualized Occlusion Combos are offered in two combinations and available in 12 popular Vita type shades. ViForm teeth were created for the professional who prefers Vita shades and Bio-


form / Portrait molds. All of our Justi Hollywood Teeth brands are crafted at American Tooth Industries's Oxnard, California factory.



JUSTI BLEND® *Lingualized Occlusion Combos™*

Justi Blend lingualized occlusion combos are offered in six combinations. The use of lingualized occlusion offers the prospect of greater patient comfort and chewing efficiency. We recommend that the mounting and ar-

ticulation of these combinations be done by specialists. See the Justi Blend mold chart for other molds.





Justi Products
A Brand Division of American Tooth Industries
 1200 Stellar Drive • Oxnard, CA 93033-2404
 (805) 487-9868 • Toll Free (800) 235-4639 • FAX: (805) 483-8482
 E-mail: info@americantooth.com • Web Site: www.americantooth.com

Order Form

Justi Lingualized Occlusion Combos™

Imperial® / ViForm® / Justi Blend® Posterior Teeth

Order # _____ Cust.# _____ Date _____
 Bill To _____ Ship To _____
 Address _____ Address _____
 City _____ State _____ Zip _____ City _____ State _____ Zip _____
 Phone # _____ Fax # _____ Phone # _____ Fax # _____
 Email _____ Signature _____

Imperial Combinations 2X8			Imperial Shades																		Imperial Teen and Young Adult Shades			
Product Number	Upper 33° Molds	Lower 10° Molds	20	25	35	45	50	59	60	62	64	65	66	67	68	69	70	77	81	82	17	18	19	21
I COMBO 1	32 S	32 X																						
I COMBO 2	30 M	30 X																						
I COMBO 3	30 S	30 X																						
I COMBO 4	34 S	34 X																						
I COMBO 5	34 M	34 X																						
I COMBO 6	32 L	32 X																						
I COMBO 7	30 L	30 X																						
I COMBO 8	32 M	32 X																						
I COMBO 9	34 L	34 X																						
TOTAL																								



ViForm Combinations 2X8			ViForm/Vita Shades											
Product Number	Upper 33° Molds	Lower 10° Molds	59 A1	62 A2	65 A3	70 A3.5	81 A4	20 B1	66 B2	77 B4	20 C1	20 C2	C3	C4
V COMBO 1	30 M	30 F												
V COMBO 2	32 M	32 F												
TOTAL														



Justi Blend Combinations 2X8			Justi Blend / BioBlend Shades								
Product Number	Upper 33° Molds	Lower 10° Molds	200 100	202 102	204 104	206 106	208 108	210 110	212 112	214 114	216 116
JB COMBO 1	32 L	32 F									
JB COMBO 2	34 L	34 F									
JB COMBO 3	32 M	32 F									
JB COMBO 4	34 M	34 F									
JB COMBO 5	30 M	30 F									
JB COMBO 6	30 L	30 F									
TOTAL											



Imperial 2x8 Posterior Tooth Cards

Number of Combos: _____ x2

1 X 8 CARDS
 GRAND TOTAL: _____

Item# 083-50-020

ViForm 2x8 Posterior Tooth Cards

Number of Combos: _____ x2

1 X 8 CARDS
 GRAND TOTAL: _____

Item# 083-50-010

Justi Blend 2x8 Posterior Tooth Cards

Number of Combos: _____ x2

1 X 8 CARDS
 GRAND TOTAL: _____

Item# 079-58-0003

1x16 or 2X8 Set - Occlusion Combination Posterior Teeth in Premium Packaging (The Combo)

Selected posterior teeth having different cusps angles to achieve the Lingual Occlusal effect for superior mastication. Our expertly crafted hard copolymer teeth are cross-linked and color stable. These selected posterior combinations are available in the full range of shades.

BioBlend® and Vita® are not a trademark of American Tooth Industries

