



Since 1864

ACRYLICS

Justi Denture Base

Microlink Powder & Liquid

High Impact for Traditional and Miracleflask Microwave Techniques

A high grade methacrylate polymer combined with a cross-linked, solvent-resistant monomer, has advantages not found in ordinary denture base materials. It has a short setting and long working time. It is hard, dense and structurally strong, having good model adaptation and high color fidelity. The high quality cross-linked monomer prevents crazing from solvents such as ether, alcohol, monomer, acetone, etc. This material will bond with any of Justi's tooth lines.

GROUP/ITEM # DESCRIPTION

Justi Denture Base Kit (20 units)*

- 1 lb. powder and 8 oz. clear liquid
- 072-04-0001 Regular pink powder
 - 072-04-0006 Dark powder
 - 072-04-0009 B Shade powder
 - 072-04-0012 Regular fibered powder
 - 072-04-0015 Dark fibered powder
 - 072-04-0018 Clear powder
 - 072-04-0033 Ethnic Natural
 - 072-04-0037 Ethnic "A"
 - 072-04-0041 Ethnic "B"
 - 072-04-0045 Regular Original Fibered

Refills - 1 lb. powder

- 072-04-0002 Regular pink powder
- 072-04-0007 Dark powder
- 072-04-0010 B Shade powder
- 072-04-0013 Regular fibered powder
- 072-04-0016 Dark fibered powder
- 072-04-0019 Clear powder
- 072-04-0034 Ethnic Natural
- 072-04-0038 Ethnic "A"
- 072-04-0042 Ethnic "B"
- 072-04-0046 Regular Original Fibered

Refills - Liquid

- 072-04-0003 Denture Base Microlink Liquid 8 oz.
- 072-04-0031 Denture Base Microlink Liquid 1 qt.
- 072-04-0005 Denture Base Microlink Liquid 1 gal.
- 072-04-0021 Denture Base Micro-Liquid 8 oz.**
- 072-04-0032 Denture Base Micro-Liquid 1 Qt.**
- 072-04-0022 Denture Base Micro-Liquid 1 gal.**

GROUP/ITEM # DESCRIPTION

Refills - 7 lb. bulk powder

- 072-04-0025 Regular pink powder
- 072-04-0026 Dark powder
- 072-04-0027 B Shade powder
- 072-04-0028 Regular fibered powder
- 072-04-0029 Dark fibered powder
- 072-04-0030 Clear powder
- 072-04-0035 Ethnic Natural
- 072-04-0039 Ethnic "A"
- 072-04-0043 Ethnic "B"
- 072-04-0047 Regular Original Fibered

Refills - 25 lb. bulk powder

- 072-04-0004 Regular pink powder
- 072-04-0008 Dark powder
- 072-04-0011 B Shade powder
- 072-04-0014 Regular fibered powder
- 072-04-0017 Dark fibered powder
- 072-04-0020 Clear powder
- 072-04-0036 Ethnic Natural
- 072-04-0040 Ethnic "A"
- 072-04-0044 Ethnic "B"
- 072-04-0048 Regular Original Fibered

Shade Guide 10 Tabs

- 072-04-0102 Justi Denture Base Shade Guide - 10 Shades

All Justi denture base products can be used with microwave curing techniques.

*Enough material for 20 dentures **For faster 3-minute curing use Justi Micro-liquid monomer. ***Customized Shade Guide available

Major Fibers



Synthetic fibers for total and partial prosthesis.

Group/Item# 073-C9000 Description Fibers 5-grams

New Shade Available 'Regular Original Fibered'



***Shade Guide Available Great Marketing Tool!

Justi Repair Material

Available in 3 shades of powder: Regular Pink, Regular Fibered Pink and Clear.

A universal type repair material, applicable for repairing all types of acrylic dentures, temporary prosthesis and indirect relines. This self-curing methyl-methacrylate material is cross-linked and internally pigmented to give a dense, even texture.

GROUP/ITEM # DESCRIPTION

Justi Repair Material Kit (20 units)

- 1 lb. powder and 8 oz. liquid
- 072-06-0001 Regular pink
 - 072-06-0004 Regular fibered pink
 - 072-06-0006 Clear

Shade Guide

- 072-06-0100 Repair Material Shade Guide

Refills - 1 lb. powder

- 072-06-0002 Regular pink
- 072-06-0005 Regular fibered pink
- 072-06-0007 Clear

Refills - 7 lb. bulk powder

- 072-06-0014 Regular pink
- 072-06-0016 Regular fibered pink
- 072-06-0017 Clear

Refills - 25 lb. bulk powder

- 072-06-0010 Regular fibered pink
- 072-06-0011 Regular pink
- 072-06-0012 Clear

Refills - Liquid

- 072-06-0003 Liquid 8 oz.
- 072-06-0015 Liquid 1 qt.
- 072-06-0013 Liquid 1 gal.



Justi Study Model Material

A methyl methacrylate self-curing material that can be poured into a silicone or rubber mold with a recommended technique that will duplicate the model in an off-white/cream color. It can be used for fabricating study model cases that emphasize the appliance rather than the model. This material can be used for many other applications where an off-white resin is needed with an autopolymerization system. Cure time is approximately 20 minutes in open air. Pressure pot technique is recommended for a more compact, dense model.

GROUP/ITEM

- 072-13-0001 Study Model Kit (1 lb.)
1 lb. powder (450gr.)
8 oz. Liquid (240ml.)
- 072-13-0002 Study Model Kit (5 lb.)
5 lb. powder (2250 gr.)
1 qt. Liquid (960ml.)

GROUP/ITEM

- 072-13-0003 Study Model Powder
1 lb. powder (450gr.)
- 072-13-0004 Study Model Liquid
8 oz. liquid (240ml.)



tcs[®] Removable Partial Material

Distributed by ATI



Beauty and Comfort...

Thermoplastic Comfort Systems (tcs[®]) Unbreakable Partial

For patients that do not want the unreliability of temporary appliances such as flippers, revealing partials such as acrylic with metal clasps, or expensive options like implants that require highly invasive surgery and long healing periods, tcs[®] Unbreakable partials may be the perfect choice. The tcs[®] Unbreakable partial allows the dentist to provide a lightweight flexible partial denture of exceptional comfort, beauty and reliability. tcs[®] Unbreakable is an extremely resilient high-performance material; tcs[®] restorations adjust and flex with the natural contours of the mouth.


tcs[®] Unbreakable's lightweight and very durable nylon-based material means that thinner and more comfortable dentures can be created. Patients will not experience that awkward bulkiness that makes most denture wearers feel uncomfortable with their appearance. As the name implies, tcs[®] is unbreakable. There is no need to worry if the denture falls in the sink or even on a hard floor. Yet, incredibly, it is much thinner than conventional appliances.

tcs[®] Unbreakable benefits:


- The ideal degree of flexibility.
- Non-invasive treatment.
- Biocompatible and ideal for patients that are allergic to acrylic monomer.
- Low water absorption, that means the partial will not warp nor become brittle.
- Can be used in combination with metal framework or precision attachments.
- Natural hue allows the appliance to simulate real gums and teeth more organically.
- Unsightly metal clasps are simply gone forever. tcs[®] Unbreakable uses tissue colored clasps, designed to blend with natural dentition.
- Lifetime warranty against breakage
- Made in the USA



(S) Small appliance 5 per Pkg - 7.9g/tube



(M) Medium appliance 5 per Pkg - 17.5g/tube



(L) Large appliance 5 per Pkg - 25.5g/tube

GROUP/ITEM #	DESCRIPTION	Color:
099-4002-01/S1	Removable partial material (S) shade 1	Light Pink
099-4002-01/S2	Removable partial material (S) shade 2	Standard Pink
099-4002-01/S3	Removable partial material (S) shade 3	Light/Dark Pink
099-4002-01/S4	Removable partial material (S) shade 4	Dark Pink
099-4002-01/SN	Removable partial material (S) natural	Natural
099-4000-01/M1	Removable partial material (M) shade 1	Light Pink
099-4000-01/M2	Removable partial material (M) shade 2	Standard Pink
099-4000-01/M3	Removable partial material (M) shade 3	Light/Dark Pink
099-4000-01/M4	Removable partial material (M) shade 4	Dark Pink
099-4000-01/MN	Removable partial material (M) natural	Natural
099-4005-01/L1	Removable partial material (L) shade 1	Light Pink
099-4005-01/L2	Removable partial material (L) shade 2	Standard Pink
099-4005-01/L3	Removable partial material (L) shade 3	Light/Dark Pink
099-4005-01/L4	Removable partial material (L) shade 4	Dark Pink
099-4005-01/LN	Removable partial material (L) natural	Natural

099-8015-01 tcs shade guide - 4 shades



- 1-Light Pink
- 2-Standard Pink
- 3-Light/Dark Pink
- 4-Dark Pink

FREE Shade Guide with purchase



Justi Crown & Bridge Formula



An outstanding cross-linked polymethyl methacrylate resin to formulate Crowns, Bridges, & Inlays. It is a heat cure method for permanent cases and long term provisionals which is available in 17 different dentin shades: 20, 25, 30, 35, 40, 45, 50, 55, 60, 62, 65, 66, 67, 69, 70, 77 and 81. Also in two enamels: light incisal and grey incisal and a finishing transparent shade. Complete technique instructions available on request.

GROUP/ITEM #	DESCRIPTION
070-40-0601	1 oz. "SR" Monomer (Solvent Resistant)
070-40-0604	4 oz. "SR" Monomer (Solvent Resistant)
070-40-0632	32 oz. "SR" Monomer (Solvent Resistant)
070-40-0699	1 gal. "SR" Monomer (Solvent Resistant)
070-40-0401/SHADE	1 oz. Powder - all shades

GROUP/ITEM #	DESCRIPTION
070-40-0404/SHADE	4 oz. Powder - all shades
070-40-0420/SHADE	1 lb. Powder - all shades
070-42-0421	Crown and Bridge Formula (SR Shade Selector)
070-40-0401/T	1 oz. Shade Transparent
070-40-0404/T	4 oz. Shade Transparent
070-40-0420/T	1 lb. Shade Transparent



Justi Imperial Namilon

An internally pigmented, cross-linked methyl methacrylate resin for permanent crown and bridge restorations. It is ideal for quality, high strength, long term provisional cases such as implants and full mouth reconstructions. Namilon is to be used with the crown and bridge investment procedures, i.e. wax-up, investment, boil-out and packing. Because Namilon is internally pigmented, it will not streak, bleed or change color. Namilon features an extended working time, and once polymerized in heat it cures quickly. Restorations are durable, solvent and wear-resistant.

GROUP/ITEM #	DESCRIPTION
070-20-1000	Namilon Kit Contents: 16 - 1 oz. bottles Namilon Powder in Imperial Shades: 60, 62, 64, 65, 66, 67, 68, 69, 70, 77, 81, 82, Incisal A, Incisal B, Opaque and Transparent 2 - 1 oz. bottles Liquid Resin, 2 - Mixing Jars, 1 - Measuring Spoon, 1 - Glass Pipette, 1 - Imperial Shade Guide, 1 - Namilon Instructions

070-20-0401/SHADE	1 oz. Namilon Powder- all shades
070-20-0404/SHADE	4 oz. Namilon Powder- all shades
070-20-0420/SHADE	16 oz. Namilon Powder - all shades
070-20-0601	1 oz. Namilon Liquid Resin
070-20-0604	4 oz. Namilon Liquid Resin
070-20-0608	8 oz. Namilon Liquid Resin
070-20-0632	32 oz. Namilon Liquid Resin

Justi Resin Modifying Stains

Justi Resin Modifying Stains are primarily used for modifying Justi Namilon shade powders and Justi Crown and Bridge formula. However, they can be used with other types of methacrylates. Also used for characterization of the cervical areas, for decalcification marks and check-lines in plastic teeth, etc. Available in shades: yellow, gray, white, brown, pink, blue, light orange and dark orange.

GROUP/ITEM #	DESCRIPTION
070-60-1000	Justi Modifying Stain Kit - includes 1/4 oz. of each shade
070-60-0400/SHADE	Refills - 1/4 oz. (Please specify shade when ordering)



Justi Cold-Cure Ortho Resin

An autopolymerizing material primarily used in the construction of orthodontic appliances. The powder is clear and, with the clear liquid, will produce a clear base material. The powder with the crystal liquid, will produce a crystal clear base material. A pink liquid is available. The pink liquid, with the clear powder, will produce a pink denture base type of material. The powder is a low-flow polymer that absorbs the monomer rapidly. This prevents the material from slumping when using the liquid dropper technique.

GROUP/ITEM #	DESCRIPTION	Liquid Refills	GROUP/ITEM #	DESCRIPTION	Powder Refills:	GROUP/ITEM #	DESCRIPTION
072-05-0001	Ortho Resin Kit w/Pink liquid	072-05-0003	Pink - 8 oz.	072-05-0002	1 lb.		
072-05-0006	Ortho Resin Kit w/Clear Liquid	072-05-0005	Pink - 1 gal.	072-05-0004	25 lb.		
072-05-0014	Ortho Resin Kit w/Crystal Liquid	072-05-0007	Clear - 8 oz.				
		072-05-0008	Clear - 1 gal.				
		072-05-0009	Pink - 1 quart				
		072-05-0010	Clear - 1 quart				
		072-05-0008	Crystal - 8 oz.				
		072-05-0009	Crystal - 1 quart				
		072-05-0010	Crystal - 1 gal.				



Justi Tray Material



A self-curing type acrylic for construction of precision custom trays and stabilized Occlusal Bite Rims. Justi Tray material is a coarse type powder that has a high degree of elasticity with a moderate curing time to allow ample working time. After the initial dough, the material can be manipulated without sticking to the hands. Having a dead pack, it can be adapted to the model without any significant shrinkage after the cure. The surface roughness aids in the retention of the impression material.

GROUP/ITEM #	DESCRIPTION
072-03-1030	30-Unit Professional Package
072-03-1025	Bulk Shipment
	Contents: Powder - 1-25 lb. drum
	Liquid - 1 gallon
Refills:	
072-03-0499	Powder - 1-25 lb. bulk
072-03-0699	Liquid - 1 gallon



Justi Frame Base Material

Justi Frame Base Material is a specially formulated microwave-type denture base material used in processing metal framework partial dentures. Can be rapidly cured in a microwave oven. This material can withstand higher heat curing than conventional acrylics. Used only for microwave curing with the Justi Miracleflask System. Available in two shades: regular fibered and dark fibered.

GROUP/ITEM #	DESCRIPTION
072-07-0001	Frame Base Regular Fibered - 1 lb. Kit - (1 lb. powder, 8 oz. Micro-liquid)
072-07-0002	Frame Base Dark Fibered - 1 lb. Kit - (1 lb. powder, 8 oz. Micro-liquid)



Justi Pakto Pure Alloy

Justi Pakto Pure Alloy is superior to other alloys including Vitallium. It is a non-precious alloy, cast in a perfect cylinder marked with the Justi logo. It has no porosity on the surface of the block. As it can be used in many different applications such as the frame of partials, customized frame structure for bridges and for crowns and is available at a very affordable price. The alloy disbursts insignificant amount of residues when melted, making it an ideal choice for all laboratory foundry. The Pakto alloy does not contain nickel, beryllium, or iron. For more information, please request an instruction sheet.

GROUP/ITEM #	DESCRIPTION
070-29-0001	Justi Pakto Pure Alloy (1 oz.)
070-29-0002	Justi Pakto Pure Alloy (3 oz.)
070-29-0001K	Justi Pakto Pure Alloy (1 Kilo)



Justi Denture Reliner

Semi Rigid, Soft and Very Soft

Heat Cure

An ethyl methacrylate-type material that can be heat cured traditionally or processed in the microwave oven within minutes when using the Justi Miracleflask system. This pink soft liner blends well with most denture base materials. A choice of three textures - semi rigid, soft and very soft. The material is tissue tolerable and will improve mouth comfort from spiny ridges or other causes of chronic soreness. Doctors, technicians, and patients will all be exceptionally pleased with this high quality, long lasting material.

GROUP/ITEM #	DESCRIPTION
072-14-0001	Justi Denture Reliner Kit - Semi-rigid
072-14-0002	Justi Denture Reliner Kit - Soft
072-14-0003	Justi Denture Reliner Kit - Very soft
072-14-0004	Justi Denture Reliner - 8 oz. Liquid - Semi-rigid
072-14-0005	Justi Denture Reliner - 8 oz. Liquid - Soft
072-14-0006	Justi Denture Reliner - 8 oz. Liquid - Very Soft
072-14-0007	Justi Denture Reliner - 1 lb. Powder



“The Ceramist’s Friend”

Same technique as is used with porcelains

The Justi Pakto® System is a high quality resin-to-metal alternative to Porcelain

Justi-Pakto is a new crown and bridge resin material system designed to save time and labor while offering maximum physical properties to simulate the hardness of natural dentition. Justi-Pakto is not a methyl methacrylate, a composite or a porcelain. It is 100% poly glycol dimethacrylate resin and is processed in a dry heat pressurized chamber.

The areas of application are:

- A) Crown and bridge veneers
- B) Facings on partial dentures and attachment frames
- C) Inlays of all types
- D) Maryland-type bridges
- E) Temporary crowns and bridges

Justi-Pakto can be used with any ADA recommended metal with no color change or fractures. The opaque has superior masking qualities and chemically bonds to the tooth shaded resin material.

Color verification can be done before the resin is cured since the shade is identical when cured. Justi-Pakto is shaded to match the Vita-Lumin® Shade System. (Vita-Lumin is not a registered trademark of American Tooth Industries.)

The same tools used to process porcelain can be used for Justi-Pakto. The Justi-Pakto Curing System utilizes 2 pieces of equipment:

1. A 4 quart, dry heat pressurized pot.
2. A Justi Cadco Cast Iron Range*.

* The Justi Cadco Cast Iron Range has been tested with the pressure pot system and is the recommended Pakto Pot Stove



Group/Item # 070-28-0001 - Justi-Pakto Kit

Optional Accessories



Group/Item # 070-28-0020 Pressure Pot - 4 quarts

Group/Item # 080-28-0033 Justi Cadco Cast Iron Range - Black

PRESSURE POT SYSTEM

GROUP/ITEM #	DESCRIPTION
070-28-0020	CONTENTS: 1 - Pressure Pot - 4 quarts, 1 - Curing Vessel, 2 - Security Spring Clips, 1 - Air Supply Hose (3 month warranty on system)

PRESSURE POT SYSTEM STOVE

GROUP/ITEM #	DESCRIPTION
080-28-0033	Justi Cadco Cast Iron Range - Black

JUSTI PAKTO KIT

GROUP/ITEM #	DESCRIPTION
070-28-0001	16 - 1 oz. Bottles Justi Pakto Powder in Vita® Shades: A1, A2, A3, A4, B2, B3, B4, C1, C2, C3, C4, D2, D3, Body Shades, Transparent Overlay and 2 Bottles Incisal
8 - 1/4 oz. Bottles Justi Pakto Modifying Stains: white, brown, yellow, gray, blue, pink, orange and dark orange	
2 - 1/4 oz. Bottles Justi Opaquer Powder (preshaded)	
2 - 1 oz. Bottles of Justi Pakto Liquid, 1 - 1/4 oz. Retention Beads, 1 - 1/4 oz. Bead Adhesive,	
2 - Dappen Dishes, 1 - Glass Pipette, 1 - Plastic Scoop, 1 - Instruction Manual. (Dappen dishes, eye dropper and plastic scoop are not available for re-order)	
1 - Justi Pakto Shade Guide - 13 shades, double-sided blue frame	

JUSTI PAKTO REFILLS

GROUP/ITEM #	DESCRIPTION
070-28-0002/SHADE	1 oz. Pakto Powders - All Shades (Includes Body Shades, Incisal or Transparent)
070-28-0003/SHADE	4 oz. Pakto Powders - All Shades (Includes Body Shades, Incisal or Transparent)
070-28-0004/SHADE	1/4 oz. Pakto Modifying Stains - All Shades
070-28-0005	1/4 oz. Pakto Opaquer - Preshaded

PRESSURE POT SYSTEM REPLACEMENT PARTS

GROUP/ITEM #	DESCRIPTION
072-28-0021	Curing Vessel Pot
072-28-0022	Security Spring Clip
072-28-0023	Air Supply Hose
800-T1325-00	Electric Heater w/Thermostat - 220 Volts - for Europe and export (not pictured)
072-28-0026	Gasket for Pakto Pot
072-28-0027	Replacement Handle for Pakto Pot - Top
072-28-0027/B	Replacement Handle for Pakto Pot - Bottom
072-28-0031	220V. European Voltage Converter for Cast Iron Range

JUSTI PAKTO SHADE GUIDE

GROUP/ITEM #	DESCRIPTION
070-28-0030	Justi Pakto Shade Guide - 13 shades on key ring
070-28-0032	Justi Pakto Shade Guide - 13 shades, double-sided blue frame

JUSTI PAKTO SUPPLIES

GROUP/ITEM #	DESCRIPTION	GROUP/ITEM #	DESCRIPTION
	Pakto Liquid		Justi Pakto Pure Alloy
070-28-0006	1 oz.	070-29-0001	1 oz.
070-28-0007	4 oz.		Justi Pakto Pure Alloy
	Retention Beads	070-29-0002	3 oz.
070-28-0008	1/4 oz.		Justi Pakto Pure Alloy
	Bead Adhesive	070-29-0001K	1 Kilo
070-28-0009	1/4 oz.		Pakto Technique video
		077-66-0003	

Veneering Material

Pyroplast, a unique veneering material that has unmatched physical and aesthetic properties, represents the closest step yet toward recreating the vitality of natural dentition. With Pyroplast, you will enjoy the benefits of fabricating single or multiple units in crown and bridge work, with ease and minimal time, while achieving a desirable restoration. Pyroplast material is a polyglycol dimethacrylate that provides superior strength and dental aesthetics.

Resin Veneering - Excellent physical and aesthetic properties. A build-up type material with no wax-up, no investing and no boil-out. Available in shades 59, 60, 62, 64, 65, 66, 67, 68, 69, 70, 77, 81, 82, Incisal A, Incisal B and Transparent Overlay.

Pyropaquer - A paint-on opaquer for masking metal in veneer crowns. In shades white, brown, yellow, gray, pink and preshaded (shade 62).

Pyrostain - For personalizing veneer jackets. In shades white, brown, yellow, gray, blue, pink, light orange and dark orange.

Pyroplast Technique - Pyroplast has antiquated all routine handling procedures for standard plastic veneering materials. With special accurately controlled low heat furnace (Pyro-Oven) and a temperature controlled forced air blower (Pyro-Heater), the veneer build-up is similar to porcelain. The pyroplast technic eliminates the time consuming procedures of wax-up and flasking. The veneering material is applied directly to metal and cured by accessory equipment, producing a restoration of optimum quality in minimum time.



Pyroplast kit

JUSTI PYROPLAST KIT

GRP/ITEM #	DESCRIPTION
070-10-1000	Pyroplast Kit in Imperial Shades Contents: 16 - 1 oz. bottles Pyroplast Powder in Imperial Shades: 59, 60, 62, 64, 65, 66, 67, 68, 69, 70, 77, 81, 82, Incisal A, Incisal B and Transparent overlay. 6 - 1/4 oz. bottles Pyropaquer in shades white, brown, yellow, gray, pink, and pre-shaded (shade 62). 8 - 1/4 oz. bottles Pyrostain in shades white, brown, yellow, gray, blue, pink, light orange and dark orange. 2 - 1 oz. bottles Pyroliquid 1 - Imperial Shade Guide, 2 - Sable brushes (#0 and #2), 4 - Dappen Dishes with Tray, 1 - Eye Dropper, 1 - Instruction Manual The Dappen Dishes and the Sable Brushes are not available for reorder.



Pyro-Oven

PYRO-HEATER AND PYRO-OVEN

GRP/ITEM #	DESCRIPTION
080-15-0810	Pyro-Heater - 110 volts, 60 cycles (3 month warranty)
080-15-0820	Pyro-Oven - 110 volts, 60 cycles (3 month warranty)



Pyro-Heater

JUSTI PYRO-PLATE

The Justi Pyro-Plate is designed for use with the Pyroplast technique. With the development of Pyroplast Crystal Monomer and the Pyro-Plate, the Pyroplast technique is simplified even further. The Pyro-Plate becomes an essential part of the technique, because it has been designed to hold exact preheating temperatures which are vital to ensure optimum properties on a consistent basis.

GRP/ITEM #	DESCRIPTION
070-08-0001	Pyro-Plate - 120 volts, 60 cycles, 65 watts. Complete with stainless tray. (3 month warranty)
070-08-0002	Stainless steel tray
070-09-0001	PyroPlast Crystal Monomer 1/2 oz.



Pyro-Plate

JUSTI PYROPLAST REFILLS AND BULK ORDERS

GRP/ITEM #	DESCRIPTION
070-10-0401/SHADE	Pyroplast Powder - all shades 1 oz.
070-10-0404/SHADE	Pyroplast Powder - all shades 4 oz.
070-11-0400/SHADE	Pyropaquer - white, brown, yellow, gray, pink and preshaded (shade 62) - 1/4 oz.
070-11-0401/SHADE	Pyropaquer - 1 oz.
070-12-0400/SHADE	Pyrostain - white, brown, yellow, gray, blue, pink, light orange and dark orange - 1/4 oz.
070-10-0601	Pyroliquid - 1 oz.
070-10-0604	Pyroliquid - 4 oz.
070-10-0608	Pyroliquid - 8 oz.

Justi Opaquer



Since 1864

A resin opaquer for masking metal in veneer crown with outstanding quality of adhesion. It is also very helpful in holding veneer color. Primarily used with Namilon or Justi Crown and Bridge Formula. However, it can also be used with other polymer-based crown and bridge materials. Brown and pink opaquers are used when experience shows they are needed.

GROUP/ITEM #	DESCRIPTION
070-30-1000	Justi Opaquer Kit Contents: 2 - 1/2 oz. liquid 5 - 1/4 oz. Powder in shades white, gray, yellow, brown and pink 1 - Sable application brush 1 - Glass pipette 1 - Polyethylene mixing slab Instructions for use
070-30-0400/SHADE	Powder Refills - all shades - 1/4 oz.
070-30-0600	Liquid Refills - 1/2 oz.



Justi Resin Cement

For Permanent Cementation

A unique copolymer offering greater adhesive properties and insolubility. Justi Resin Cement offers a lower film thickness and is radio-opaque.

ADVANTAGES:

- Shows adhesion to tooth structure superior to any other cement tested after 24 hours immersion in water.
- Meets the requirement for film thickness in ADA Specification #8, when used according to instructions.
- Available in six shades - Natural Yellow, Natural Gray, White, Pink, Brown and Clear.
- Insoluble in mouth fluids.
- As harmless to the pulp as present resin filling materials.
- Radiopaque.

GROUP/ITEM #	DESCRIPTION
070-80-1000/SHADE	Justi Resin Cement Kit Contents: 2 - 1/2 oz. Powder; 1 - 1/2 oz. Liquid; a Measuring Spoon, 1 Eye Dropper and Instructions for Use
070-80-0400/SHADE	Powder Refill (1/2 oz.) and a Measuring Spoon.
070-80-0600	Please specify shade when ordering Liquid Refill (1/2 oz.) and Eye Dropper



Justi Temporary Crown & Bridge Resin - Chairside Product

3 Shades - Light, Medium and Dark

An autopolymerizing self-cure powder and liquid resin for fabricating a crown, bridge or splint. It is very hard, color stable, easy to polish and comes in light (shade 59), medium (shade 62), and dark (shade 69). This material has a low exothermic heat peak, and therefore, is tolerable to the oral cavity; it does not burn.

GROUP/ITEM #	DESCRIPTION
072-25-0001	Professional Kit - Boxed 3 oz. each light, medium and dark powders and 8 oz. liquid
Refills:	
072-25-0002	Light Powder (85 grams) 3 oz.
072-25-0003	Medium Powder (85 grams) 3 oz.
072-25-0004	Dark Powder (85 grams) 3 oz.
072-25-0005	Clear Liquid (240 ml) 8 oz.



Justi Tooth Shades - Laboratory Product

Cold Cure Temporary Crown and Bridge Acrylic for Labs

Specially prepared formulation for dental laboratories to construct and build strong temporary crowns and bridges. Crosslinking monomer with high craze resistance and toughness. Kit available in Imperial/Bioform type shades: 59, 62, 65, 66, 67, 69, 81. **New Teen Shades: 17, 18, 19, and 21 available in 1lb. bottle only.**

GROUP/ITEM #	DESCRIPTION	Refills available:
072-25-0010	Complete Tooth Shade Kit Contents: 7 -3 oz. tooth shade powders, 2 - 8 oz. bottle of monomer. 2 - Dappen dishes and eye dropper	072-25-0008 8 oz. bottle of monomer 072-25-0009/SHADE 3 oz. bottle of tooth shade powder Include shade number with request
072-25-0011/shade	Imperial Teen Shade Powder Contents: 1 Lb. bottle of tooth shade powder, be sure to indicate the shade - example 072-25-0011/18.	

Note: Keep material away from heat and light to maintain a long shelf life.



Resin Instrument Kit



Instrument set for resins and composites modeling.

GROUP/ITEM #	DESCRIPTION
073-B9000	Resin Instrument Kit



Light Cure Trays

Light-curing plates for modeling individual impression trays and base-plates. Available in standard pink and transparent, in the upper form only.

GROUP/ITEM #	DESCRIPTION
073-C1501	Light Cure Tray (50 pcs.) Pink
073-C1502	Light Cure Tray (50 pcs.) Transparent



Itoya Pen

The Itoya Permanent Black Ink Pen is the only pen that will easily write on the most polished acrylic materials. Non smearing. Available for the first time as a separate item.

GROUP/ITEM#	DESCRIPTION
077-65-0054	.1 Itoya Permanent Black Ink Pen



STG System



This unique system provides detailed characterization of the gum tissue in prosthodontic appliances. A must have for the consummate technician.

GROUP/ITEM#	DESCRIPTION
071-67-0001	STG System: Resin Stains for Artificial Gums
071-67-0002	Gums Shade Guide Plus Metal Mold to Make Gums
071-67-0003	PALATO – Denture Base
071-67-0006	Instructional Video (Portuguese) (Written instructions in English included)

