

ASA DENTAL PRODUCTS

DISTRIBUTED BY AMERICAN TOOTH INDUSTRIES

Facebow and Transfer Table

The “Easy Bow” Facebow transfer system

Because Every Patient is Different...

The “Easy Bow” 5020 transfer system set includes the facebow with three forks and a transfer table to move patient articulation characteristics to an ASA 5032 articulator. An instructional DVD and printed instructions in English, French, Spanish and Italian are included in the kit. Recommended for Prosthodontists and Dental Laboratory Specialist. Made by ASA for American Tooth Industries.



GROUP/ITEM # DESCRIPTION

098-5020 Facebow & Transfer Table

High Tech 5032 Articulator

A lightweight anatomic articulator made of die-cast aluminum. This quality instrument has many ease-of-use design features for a quick and accurate model fixation. The 5032 high-Tech articulator is designed to work seamlessly with the 5020 “Easy Bow” facebow transfer system.

- Adjustable back rod by registration screw without incisal rod by keeping constant trajectories and inclinations during side and back/forth movements.
- Automatic return to hinge axis with springs.
- Anti-shock hinge rod in case articulator accidentally falls.
- Big size lock screws for mounting plates.

English instructions brochure is included in box. Made by ASA for American Tooth Industries.

GROUP/ITEM # DESCRIPTION

098-5032 High Tech Articulator

