




Since 1864

DENTURE ACCESSORIES

Denture Specialist Reception Brochure

A great idea for the dental reception area. This brochure will inform new denture wearers about the proper care and maintenance of their new dentures. Available plain or imprinted with personalized information.



GROUP/ITEM #	DESCRIPTION
077-65-0044	Denture Specialist Reception Brochure - not imprinted
077-65-0044/P	Denture Specialist Reception Brochure w/custom imprint

Certificate of Authenticity / Warranty Card


Prosthetic Appliance; Quality Assurance Inserts for Justi tooth lines. These handy Certificates, which now includes the Justi Warranty Card, are available for you to prove to your patients that the finest quality of teeth have been used in the fabrication of their prosthetic appliance.



GROUP/ITEM #	DESCRIPTION
077-65-0021	Certificate of Authenticity/ Warranty Card 20 sheets per pack

Justi Color Neutralizer

Justi Products invented this accessory to be used as a neutral background for choosing the correct shade from Justi teeth shade guides. The neutral blue color enables the eye to focus on the shades without distraction from surrounding colors.



GROUP/ITEM #	DESCRIPTION
077-65-0037	Justi Color Neutralizer

Kroma Tabs

Justi Products invented this accessory to be used in the patient's mouth as a neutral background for choosing the correct shade from the Justi teeth shade guides. The neutral blue background enables the eye to focus on the shades, and not be distracted by surrounding colors. U.S. Patent #5,257,931.



GROUP/ITEM #	DESCRIPTION
086-51-0415	Kroma Tabs 50/Pk.

Denture Bags




Place Prosthodontic and/or Orthodontic appliances in these strong, transparent bags with your choice of disinfecting solution prior to a hermetic seal. Bags available with or without personalized imprint.

GROUP/ITEM #	DESCRIPTION
077-65-0050	1000 Denture Bags (no imprint)
077-65-0051	5000 Denture Bags (imprinted)

Rx Pads

Remind your customers about you each time they write a prescription. Customized prescription pads are a simple and effective way to increase your business.



GROUP/ITEM #	DESCRIPTION
077-65-0005	Prescription Pads (Our design/Your Logo 100 pads (20 sheets/pad)
077-65-0005	Prescription Pads (Your design/Your Logo 100 pads (20 sheets/pad)

Customized Shade Guides



Justi Imperial, Justi Dymon-Hue HPT, Justi Blend, ViForm, Regal-D-Blend and Vita-Type Teeth.

To improve shade communications, and provide an outstanding marketing tool to dental labs and clinics, Justi has created the Customized Shade Guide complete with the Justi Eterna Shield. U.S. Patent # 5,257,931. With a minimum order of 30 shade guides of one type, customized shade guides can be imprinted with your name and a marketing message at the regular shade guide price. Customized shade guides allow the option of tailoring both the front and back message areas and the number of shade tabs. The message areas are free form and can contain name, address, shade information, and marketing message limited only by imagination and the available text space. Tab selection for a customized shade guide can be from 1 to the maximum number of available shades or any number in between (Example: To order a 9 tab shade guide from a standard Justi VI-Type shade guide, which has 16 shades, on the order form simply specify the 9 shade tabs desired). The one exception is the Imperial shade guide which is only available with 18 shades.

Now you can build name recognition for your lab and make more money selling Justi teeth in your denture, overdenture and partial cases. It is recommended that the shade guide be disinfected with Cidex Plus Solution.

GROUP/ITEM #	DESCRIPTION
084-50-0125	ViForm Shade Guide in 12 shades
084-50-0903	Regal-D-Blend Maxillary Centrals Shade Guide in 6 shades
084-51-0413	Justi Imperial Maxillary Centrals Shade Guide in 12 shades
084-51-0425	Imperial Shade Guide w/Kroma Tab - 18 Age Shades
084-59-0026/C	Justi Blend Maxillary Centrals Shade Guide in 9 shades
084-19975	Dymon-Hue HPT Maxillary Centrals Shade Guide in 12 shades
084-04-0165	Vita-Type Shade Guide with 12 shades
084-04-0166	Vita-Type Double Sided Shade Guide in 16 shades

The number of shades listed for each shade guide is the maximum shades available for the specified shade guide. Call for details.

Dental Saw

Fits beautifully in the hand. Handle is made of aluminum and is light weight. Casing for blade is made of stainless steel with a curvature for easily eyeing the models. Round ball at the end of the handle is used to release the tension on the blade or replace the blade. This method gives you the same tension every time you insert the blade. Quick, easy and make a nice narrow cut in the plaster, stone or other gypsum products. For general use.

GROUP/ITEM #	DESCRIPTION
027-20027	Dental Saw - includes 1 blade
027-20028/L	Saw Blades, Large 12/pkg



Reliner


Adjustable appliance for cold rebasings in general.

GROUP/ITEM # 073-E3100



IsolMajor


Universal plaster/resin separator. 500 ml.
GROUP/ITEM # 073-L2710 IsolMajor 500 ml.



Facial Indicator

Used to determine anterior tooth size and mold in total prosthesis.

GROUP/ITEM #
074-T7900 Facial Indicator



I-Shade Guide

A unique system for the construction of individual shade guides using acrylic or composite materials.

GROUP/ITEM # DESCRIPTION
073-S6000 I-Shade Guide - 40 Tab Shade Guide to customize

