

# SOLDER & ACCESSORIES

## Soldering Tweezers

Stainless steel soldering tweezers with plastic handles for protection from heat of soldering. Available in two styles: straight or curved tips.



GROUP/ITEM #	DESCRIPTION
007-6400-1	Soldering Tweezers Straight Tip 6 1/2"
007-6400-2	Soldering Tweezers Curved Tip 6 1/4"

## Solders for Stainless Steel



For use with bridges and crowns.  
Cadmium free.  
Packs of 10 g approx.

ITEM # R0220-00  
Solder for Stainless Steel  
Plate (0.43 mm thick)

ITEM # R0221-00  
Solder for Stainless Steel  
1 mm diameter Wire

Leone has been Manufacturing Silver Solder since 1940.

	Melting Point	Silver Content	
	750° C 1382° F FRAMEWORK	27% AG	
	750° C 1382° F FRAMEWORK	27% AG	6 1/2 dwt.

## Orthodontic Solders

For all orthodontic brazing work.  
Cadmium free. Pack of 10 g approx.

ITEM # R0222-00 Wire Form, 0.40 mm dia.

ITEM # R0223-00 Tape Form, 1 x 0.5 mm.

ITEM # R0224-00 Rods, 1 mm dia.

Leone has been Manufacturing Silver Solder and Flux since 1940.

	Melting Point	Silver Content	
	600° C 1112° F ORTHO	65% AG	6 1/2 dwt.
	600° C 1112° F ORTHO	65% AG	6 1/2 dwt.
	600° C 1112° F ORTHO	65% AG	6 1/2 dwt.



## Thermo Insulator Paste for Solders

It protects plaster and resin during soldering procedure. Excellent consistency.  
Easy to remove with water.  
Pack of 235 mL

ITEM #  
R0227-01 Thermo Insulator Paste



# SOLDER & ACCESSORIES

## *Flux*

In paste, for solders.  
Pack of 30 g approx.

ITEM #  
R0226-00

Flux



## *Solder Powder*



Powder for solders. Pack of 50 g  
approx.

ITEM #  
R0225-00 Solder Powder