

JUSTI TOOTH TRAYS, CABINETS & CARRYING CASES

Justi Tooth Trays and Cabinets are an ideal Storage and Inventory System for Justi Teeth. Cabinets are available in three sizes.

- Large will hold up to 4576 combination Justi anterior and/or posterior tooth cards.
- Medium will hold 1512 anterior and/or posterior cards.
- Small has a capacity of 210 to 360 cards of anteriors and/or posteriors.

The large counter cabinet does not require inserts. Individual tooth drawer trays contain inventory strips to write the card brand, shade and mold information which simplifies reorders. Justi's tray inserts also fit a standard Dentsply 16 1/2" x 16 7/16" metal tray which minimizes brand conversion costs.

44 Drawer Tooth Counter

(Does not include inserts)


Capacity:	Tooth Capacity	Total Cards
Tooth Line	Cards per Drawer	Per Counter
Dymon-Hue HPT Anteriors	72	3,168
Dymon-Hue HPT Posteriors	104	4,576
Justi Imperial Anteriors	80	3,520
Justi Imperial Posteriors	104	4,576
Justi Blend Anteriors	64	2,816
Justi Blend Posteriors	104	4,576

Counter Dimensions: 19" x 19" x 47 1/4"*

*Please allow an additional two inches when casters are installed.

GROUP/ITEM #	DESCRIPTION
086-55-0200	44 Drawer Tooth Counter
086-55-0101/T8	Tooth Counter Tray 8 Slots

Eight sections per drawer.
Drawers are made of
hard compressed plastic.



27 Drawer Lab Tooth Counter


Laminated Material. 27 Drawers with a capacity of 1512 anterior and/or posterior cards. Can be used for any Justi tooth brand.

Counter Dimensions: 19" x 19" x 33 1/2"*

*Please allow an additional two inches when casters are installed.

GROUP/ITEM #	DESCRIPTION
086-55-0201	27 Drawer Lab Tooth Counter
086-55-0201/T	Tooth Counter Trays 10 & 4 Slots

With each drawer we supply 8 Universal Insert Trays for a total storage of 56 cards. Drawers are constructed of hard compressed plastic. Four sections per drawer.



Justi Lab Tooth Cabinet

Made of heavy plastic. Ten drawers with 3 universal insert tooth trays each. Will hold 210 to 360 anterior and/or posterior tooth cards. Can be used for all Justi tooth brands.



Cabinet Dimensions:
9 1/2" x 12 7/16" x 11 5/16"

GROUP/ITEM # DESCRIPTION
086-55-0101 Lab Tooth Cabinet



Drawer with Universal Trays Inserts


Tooth Brand Conversions are Easy!

You can use the standard Dentsply metal tray (16 1/2" x 16 7/16") with these inserts and store up to 56 cards in each drawer!

Justi Lab Tooth Cabinet Drawer Inserts

GROUP/ITEM #	DESCRIPTION	GROUP/ITEM #	DESCRIPTION
077-65-0007	Tooth Drawer Insert - holds 20 sets of the following: Justi Blend Anteriors, ViForm Anteriors and/or Posteriors, Major Plus Comp Posteriors and Triplex Posteriors.	077-65-0027/S*	Tooth Strip Insert - for 36 Set Tooth Drawer, Anteriors.
077-65-0007/S*	Tooth Strip Insert - for 20 Set Tooth Drawer, Anteriors and/or Posteriors.	077-65-0028	Tooth Drawer Thick Insert - for Posteriors. Holds 28 sets of the following: Justi Imperial Posteriors, Justi Blend Posteriors, Dymon-Hue HPT Posteriors and Justi Imperial Block Posteriors.
077-65-0015	Universal Tray Insert - fits 7 Anteriors and/or Posteriors. Fits all Justi tooth cards.	077-65-0028/S*	Tooth Strip Insert - for 28 Set Tooth Drawer, Posteriors.
077-65-0016*	Tooth Strips Universal - 7 spaces for tooth card information.		
077-65-0027	Tooth Drawer Thick Insert - for Anteriors. Holds 36 sets of the following: Justi Imperial Anteriors, Dymon-Hue HPT Anteriors, Super Lux Anteriors, Bambino and Justi Kind Anteriors and Posteriors and Major Dent Anteriors.		* Tooth Strip Inserts are sold separately

Tooth Drawer Insert
077-65-0007



Universal Tray Insert
Fits 7 Anteriors or 7 Posteriors
077-65-0015



Tooth Drawer Thick Insert
For Anterior Teeth
077-65-0027



Tooth Drawer Thick Insert
For Posterior Teeth
077-65-0028



Drawer inserts shown with optional Tooth Strips.

TOOTH TRAYS, CABINETS & CARRYING CASES

Pozzi Transport Case

This vinyl transport case with castors has been specially designed to perform as a small, roomy "take-me-along" orthodontic office. It can be easily moved without damage to the materials stored inside the compartments.

GROUP/ITEM #	DESCRIPTION
019-T0020-01	Pozzi Empty Transport Case with Trays/Inserts on Wheels

