

WAXES & SOLVENTS

Justi Wax Solvent No. 20

An efficient, safe and pleasant scented solvent used for wax elimination and shining after waxing a prosthetic case. An excellent separator for plaster, stone and any other gypsum product. Laboratories have reported to us that they also use Justi Wax Solvent #20 in the following applications:

1. Added to boil tank water.
2. To clean counter tops.
3. To clean tools and equipment.
4. To remove wax from clothing.

GROUP/ITEM #	DESCRIPTION
070-01-0616	1 Pint
070-01-0632	1 Quart
070-01-0699	1 Gallon



Pozzi Orthodontic Relief Wax

Packaged in a convenient patient packet containing five pieces of natural, protective relief wax. It is specially manufactured for orthodontic use. It protects brackets and tubes during the bonding procedure. It may also be used on brackets and ligature wires for protection of the mucous membrane against sharp or irritable objects in the mouth. The natural composition of the wax makes it absolutely safe if accidentally swallowed. *Available personalized upon request.*

GROUP/ITEM #	DESCRIPTION
066-66-004RW	Orthodontic Relief Wax - 50 packets/box



Pozzi Orthodontic Bite Wafers

These bite wafers were designed for taking accurate occlusal impressions for the construction of orthodontic appliances. They come in a one pound box of 24 red-tinted wafers.

GROUP/ITEM #	DESCRIPTION
066-66-003OWB	Orthodontic Bite Wafer 24/box



NEW Ortho & Utility Waxes

Pozzi Ortho Set-Up Wax

A strong, tough, dependable wax used in Orthodontics to show patients projected tooth movement before, during and upon completion of treatment. Ortho Set-up Wax melts easily, flows smoothly and cools quickly making the wax-up smooth and easy. The unique characteristics of this wax will allow you to warm the wax model in a water bath and move the teeth into the desired position without the worry of teeth "popping" out. Waxes available in sheets or rods.

GROUP/ITEM #	DESCRIPTION
066-66-007	Ortho Set-up Wax Sheets 1 lb. box
066-66-008	Ortho Set-up Wax Rods 24 rods to a Box



Justi Red/White Large Utility Wax

A soft, pliable, tacky wax that is strong enough at room temperature to allow a strip to be torn away without breaking yet as minimal amount of stretching. It adheres and adapts to most any material, yet it can be easily pulled away. Utility wax can be used to, line an impression tray, box an impression or used as a relief wax for orthodontic patients. Waxes available in red or white large strips.

GROUP/ITEM #	DESCRIPTION
066-66-009	Red Large Utility Wax Strips-75 strips to a Box.
066-66-010	White Large Utility Wax Strips-75 strips to a Box.



Justi Universal Set-up Wax

A tough yet dependable wax that melts easily, flows smoothly, and cools quickly. Can be used in a variety of intraoral applications; Set-ups, impressions, bite registration and as a soft barrier for appliances.

GROUP/ITEM #	DESCRIPTION
066-66-011	Universal Set-up Wax in sheets -1lb Box

