

# AMERICAN DENTAL SUPPLY PRODUCTS

DISTRIBUTED BY ATI

## SPECIALIZED DENTURE LABORATORY SUPPLIES

Regular &  
Heavy Duty

### Boil Out Gloves

#### Extreme High Temperature Protection for your Hands

New, extreme high temperature protection for your hands! Our new, Heavy Duty Boil Out Gloves protect your hands from hot flasks and liquids. Withstands hot metal up to 480°F. Excellent dexterity. Positive nonslip grip, puncture and great chemical resistant insulated, neoprene on cotton interlock with cotton fleece inner liner. Our Standard Boil Out Gloves are effective to 212°F. Durable and long lasting Knit lining provides heat insulation and comfort. Non-slip grip & chemical resistant.



**Group/Item# Description**

- 097-G222-1 Boil Out Gloves 12" Reg.
- 097-G222-2 Boil Out Gloves Extra Long
- 097-G222-3 Cotton Liner for Boil Out Gloves
- 097-G222-4 **New** HD Boil Out Gloves (Med 14")
- 097-G222-5 **New** HD Boil Out Gloves (Large 14")
- 097-G222-6 **New** HD Boil Out Gloves (X-LG 14")

### Pink Polishing Paste®

#### For Metal and Acrylic

Polish to a sparkling shine in seconds! Use on soft dental alloys and acrylics. Burnishes gold and shines up acrylic. Easy to use. Just impregnate felt or flannel buff wheel with Pink Polishing Paste and lightly buff at medium speed. Rines off with water.

**Group/Item# Description**

- 097-P320-1 Pink Polishing Paste  
1/2 oz. bottle
- 097-P320-2 Pink Polishing Paste  
2 oz. bottle
- 097-P320-4 Pink Polishing Paste  
8 oz. bottle



### Convection Soldering Table®

#### Professional Speed Cure Support for Place-It

Convection Soldering Turn Table guarantees perfect fitting restorations. Position work pieces on the mesh screen in the optimum position for soldering. Touch to turn. Solid refractory ceramic table/support screen combination permits Place-It to be hardened quickly, locks your long span framework accurately in register for soldering or welding. Custom supports for porcelain veneers are formed in record speed on the support screen. High density durable board radiates heat from torch back to the underside of the Place-It providing instantaneous hardening. Support screen can be removed for use on a hot plate. Support screen can be inserted directly into furnace for post-soldering. Heat resistant turntable provides superior convenience and soldering ease for several units simultaneously.



**Group/Item# Description**

- 097-S705-1 Convection Soldering Table Kit
- Kit Contains: Table, Support Frame, Two 3" mesh screens, Place-It Type2

### Diamond Polishing Paste

#### For Porcelain and Chrome

The king of abrasives! Save Money! Heavy concentration of diamond crystals in a non-greasy, water soluble base. Glaze adjustable areas and contacts on porcelain without refiring. polish hard partial alloys to a final super luster, and rinse clean.



**Group/Item# Description**

- 097-P322-1 Diamond Polishing - 2 grams
- 097-P322-2 Diamond Polishing - 5 grams
- 097-P322-3 Diamond Polishing - 10 grams

### Ty-Chrome Powder

#### Classic Universal Polish Medium

Polishes hard metal frameworks to a mirror finish. Polish acrylic restorations to a final high shine. Mixes with water or glycerine to form an economical polishing slurry.



**Group/Item# Description**

- T410-9 Ty-Chrome Powder  
(1lb tub)
- T410-10 Ty-Chrome Powder  
(5lb tub)

### Alginate Tray Cleaner

#### Save Time and Soak Trays Clean

Easy to use. Just add to water and soak impression trays clean. Removes stubborn alginate residues. Safe powder formula stores easily. Indispensable, efficient cleaner for plastic and metal trays. Clinically tested.



**Group/Item# Description**

- 097-T734-1 Alginate Tray  
Cleaner Powder - 2 lb
- 097-T734-7 Alginate Tray  
Cleaner Powder - 5 lb

### Turn N' Solder®

#### Economical, Rotating, Asbestos-Free Heat Shield

Turn N'Solder is an economical soldering board, designed for light duty, mounted in a stable, heat-resistant turntable. Super convenience and soldering ease. Pin your workpiece in place or use Place-It support. Use Turn N' Solder directly in front of porcelain furnace for cooling and manipulation of sagger trays.



Resistant to molten metals and most chemicals. Economical lightweight ceramic. Temperature tolerance of 3000°F. Withstands rapid heating and cools to touch in seconds.

**Group/Item# Description**

- 097-T703-1 Turn N' Solder Table  
(6"x 6"x 1/2")

# AMERICAN DENTAL SUPPLY PRODUCTS

DISTRIBUTED BY ATI

## SPECIALIZED DENTURE LABORATORY SUPPLIES

### Place-It Type2® & Type3®

**Place-It Type 2: Solder pontics without investing repair pin-holes in copings. Create Custom refractory pieces. Not tacky. Safe and asbestos free.**

Holds your assembly during soldering, welding, and firing. Will not rob heat from joints, insuring an excellent bond. Heat shield protects acrylic during soldering. Use to repair furnace parts. Use to 3200°F.



**Place-It Type 3: Improved for laminates & porcelain.**

Holds all ceramic and PFM restorations in positions on firing pegs. Custom support for veneers, inlays and onlays for direct firing on refractory pads. Tacky, ready to use. Asbestos free. No hot spots. Easy clean-up.

Group/Item#	Description
097-P702-5	Place-It Type2 Medium syringe 30cc
097-P702-9	Place-It Type3 Medium syringe 30cc

### Plastical Release Agent®

**Prevents Die Materials From Sticking to All Types of Impressions**

Specifically designed to work with synthetic Die Materials. A terrific separator for resins, silicones, epoxies, and composites. Easy spray on or paint. Non-flammable formula leaves a fine, thin film. Spray on steel instruments to create a non-stick surface for composite material.



Group/Item#	Description
097-P510-3	Plastic Release Agent-4 oz. spray
097-P510-4	Plastic Release Agent-8 oz. refill
097-P510-6	Plastic Release Agent-32 oz. refill

### Son-Eez Blue®

**General Purpose Ultrasonic Cleaner**

Powerful and economical! New, improved, extra-strength Son Eez Blue is a combination of ingredients safe for use on all non-porous surfaces (glass, metal, plastic, and rubber components). It has spectacular cleaning power, removing oils, organic debris, abrasives, grinding compounds, and general soil. High 1:32 dilution provides maximum economy. Use diluted cleaner directly in the ultrasonic tank as a coupler, or in a glass, or plastic beaker. promotes cavitation. Super in a heated bath with soft or hard water. Rinses clean, leaves no film. Pleasantly scented. Concentrated liquid and powder formulas available.



Group/Item#	Description
097-S924-6	Son-Eez Liq Blue 1qt
097-S924-8	Son-Eez Liq Blue 1gal
097-S924-9	Son-Eez Liq Blue 1 Lb

### Wax Shine®

**Glass finish on Dental Waxes**

Spreads easily on all dental waxes. Dries in minutes to a high shine without buffing. Masks fine scratches in the wax surface. The practical alternative to manual polishing of wax-ups, bite blocks, and wax try-ins. Forms a crystal clear, durable finish. Strengthens and hardens thin wax copings. Safe ready to use, cleans up with water.



Group/Item#	Description
097-W220-5	Wax Shine-pint
097-W220-6	Wax Shine-quart

### Solder Board®

**Convenient, Asbestos Free - Protect your Benches**

Economical protection for work tables or benches. Withstands quick heating and cools in seconds. Use in front of furnace to cool sagger trays. Use under bunsen burner to prevent scorching of formica or wood counters.



Group/Item#	Description
097-S701-1	Soldering Board 6"x6"x.5"
097-S701-2	Soldering Board 6"x12"x.5"
097-S701-3	Soldering Board 12"x12"x.5"

### Liquid-Foil Tin Foil Substitute®

**Superior Formulas Cleanly Separate Acrylic & Stone - Clear Color**

For labs that demand superior performance and quality. Ready to use. Adheres to stone, gels quickly, penetrates, and forms a fine tough film. Ideal viscosity results in smooth application whether painted on or poured. Requires only one coat. Film doesn't lift from edges. Smooth flow on stone and around teeth. Will not stick to teeth or metal. Rinses off cleanly with water. Two super stable formulas. Contains no alcohol to evaporate.



Shipped year round in any temperature.  
Premium X-FOIL guarantees a perfect separation

Group/Item#	Description
097-T550-6-CR	Original Liquid-Foil - 1 quart
097-T550-8-CR	Original Liquid-Foil - 1 gallon

### Pumice Deodorant & Disinfectant®

**Keeps Your Lab Smelling Fresh & Clean**

A unique blend of concentrated, EPA registered, biocidal ingredients designed to deodorize. Pumice Deodorant reduces microorganisms and fungi that cause foul odors and slime, keeping acrylic teeth or corrode pumice tray. Super concentrate: 1:32. Minty aroma.



Group/Item#	Description
097-P946-3	Pumice Deodorant - 4 oz. Concentrate
097-P946-6	Pumice Deodorant - Quart Concentrate
097-P946-32	Pumice Deodorant - Kit (Operator)
097-P946-30	Pumice Deodorant - Kit (lab)