



A Brand Division of American Tooth Industries


1200 Stellar Drive • Oxnard, CA 93033-2404
(805) 487-9868 • Toll Free (800) 235-4639 • FAX: (805) 483-8482
E-mail: info@americantooth.com • Web Site: www.americantooth.com



JUSTI-KURA CURING POT

The Justi-Kura™ Curing Pot has been designed to satisfy the need of the modern technician. It is sturdy, very practical, absolutely safe and has a large capacity of 6 quarts.

The Justi-Kura™ Curing Pot is designed for curing Pozzi Ortho Resin, Justi Ortho Resin, and Justi Repair Material. This system ensures a pourous-free cure. When powder and liquid application is complete, cover appliance with warm water in the pot, move the lid to the closed position and inject air 20 to 30 psi.



CHARACTERISTICS & ACCESSORIES

Table with 3 columns: No., Description, and Grp/Item #. Contains 7 rows of product specifications.

OTHER RELIABLE JUSTI PRODUCTS:

Hardened Plastic Teeth:
Justi Imperial®
Anterior & Posterior
Justi Blend®
Anterior & Posterior
Dymon-Hue® HPT
Anterior & Posterior
Miracleflask Curing System
For Dentures, Partials & Relines
Justi Ortho Resin

Justi Imperial Plastic Facings
Justi Wax Solvent
Justi Milling and Polishing Cream
Justi Resin Cement
Justi Opaquer
Pyroplast Veneering Material
Namilon Resin for Crown and Bridge
Justi Denture Base Material
Justi Repair Resin
and many more...

Box containing the text 'READ & PASS ALONG TO' followed by four horizontal lines and 'Then file for reference' at the bottom.

Additional questions, to place your order, or for the Dealer nearest you, please call the Tooth Express Order Desk at 800-628-1437.