

# Regal-D-Blend®



## *Mold Chart*

for  
Anterior and Posterior  
Hardened Plastic Teeth

ORTHOCHROMATIC  
MULTI-LAYERED  
ANTERIOR &  
INTERCUSPAL-OCCLUSAL DESIGNED  
JUSTI BLEND® POSTERIOR  
HARDENED PLASTIC TEETH  
SIX SHADE SYSTEM

Highly Cross Linked/Intersectional Polymer Network • Wear Resistant  
Vacuum Processed • Fluorescent/Photogenic  
No diatoric retention necessary.  
100% bondability with methyl methacrylate denture base material.  
Complies with ANSI/ADA Specification 15, ISO 3336, BS 3990, NF S91-213, DIN 13914





REGAL-D-BLEND MULTI-BLENDED ANTERIORS 27 UPPERS & 8 LOWERERS

<b>A4</b>			<b>4</b>
	10.0   8.3		

<b>D6</b>			<b>6</b>
	12.0   9.0		

<b>A7</b>			<b>7</b>
	10.0   8.8		

<b>E1</b>			<b>1</b>
	9.3   7.3		

<b>B3</b>			<b>3</b>
	9.5   8.0		

<b>E2</b>			<b>2</b>
	10.0   7.8		

<b>B8</b>			<b>8</b>
	11.3   8.8		

<b>E3</b>			<b>3</b>
	10.5   8.3		

<b>C1</b>			<b>1</b>
	10.0   7.5		

<b>E6</b>			<b>6</b>
	10.0   8.5		

<b>C2</b>			<b>2</b>
	10.0   8.0		

<b>E8</b>			<b>8</b>
	11.0   8.8		

<b>C5</b>			<b>5</b>
	11.8   8.8		

<b>F1</b>			<b>1</b>
	8.5   7.8		

<b>C6</b>			<b>6</b>
	10.8   8.5		

<b>F4</b>			<b>4</b>
	9.0   8.5		

<b>D1</b>			<b>1</b>
	10.0   7.5		

<b>G5</b>			<b>5</b>
	10.5   8.8		

<b>D2</b>			<b>2</b>
	10.0   8.0		

<b>G7</b>			<b>7</b>
	10.8   8.8		

<b>D3</b>			<b>3</b>
	10.5   8.3		

<b>H2</b>			<b>2</b>
	10.0   8.0		

Photographs of teeth are within plus or minus 3% of actual size.

<b>H4</b>			<b>4</b>
			51.0 44.67 — —
10.5	8.3		

<b>J2</b>			<b>2</b>
			46.5 40.75 — —
10.0	7.8		

<b>J3</b>			<b>3</b>
			48.0 42.68 — —
10.0	8.0		

<b>J5</b>			<b>5</b>
			53.0 47.14 — —
11.5	9.0		

<b>K4</b>			<b>4</b>
			49.5 43.90 — —
10.8	8.0		

<b>I</b>			<b>C1, D1 E1, F1</b>
			34.5 30.29 — —
9.0	4.8		

<b>2</b>			<b>C2, D2 E2, H2 J2</b>
			36.0 29.86 — —
8.8	4.5		

<b>3</b>			<b>B3, D3 E3, J3</b>
			37.0 30.17 — —
8.3	4.8		

<b>4</b>			<b>A4, F4 H4, K4</b>
			38.0 31.07 — —
9.0	4.8		

<b>5</b>			<b>C5, G5 J5</b>
			40.0 32.80 — —
9.8	5.3		

<b>6</b>			<b>C6, D6 E6</b>
			39.5 33.32 — —
8.5	5.3		

<b>7</b>			<b>A7, G7</b>
			41.0 33.38 — —
9.5	5.0		

<b>8</b>			<b>B8, E8</b>
			41.5 35.77 — —
9.3	5.5		

0° POSTERIORES - JUSTI BLEND 3 UPPERS & LOWERES \*

<b>29 MF</b>			6.69
			8.15
			28.5 30.5

<b>31 MF</b>			6.60
			7.62
			30.0 32.8

<b>33 MF</b>			7.31
			8.30
			32.2 34.0

10° POSTERIORES - JUSTI BLEND 3 UPPERS & LOWERES \*

<b>30 F</b>			8.41
			9.04
			30.0 33.0

<b>32 F</b>			8.72
			9.63
			31.0 33.5

<b>34 F</b>			9.15
			10.17
			33.0 36.5

20° POSTERIOBS - JUSTI BLEND 6 UPPERS & LOWERB \*




ANATOMICAL MEDIUM - 33°  
JUSTI BLEND 3 UPPERS & LOWERB \*





ANATOMICAL LONG - 33°  
3 UPPERS & LOWERB



ANTERIOR



POSTERIOR



Mold Classification Symbols: **SQUARE** (black square), **TAPERING** (black triangle), **OVOID** (black circle), **Combined Symbol** (black triangle and circle), **Tapering Ovoid Type** (black triangle and circle with a dot)

LEGEND: D - Straight line measurement of 4 upper posteriors. E - Straight line measurement of 4 lower posteriors. **JUSTI BLEND POSTERIOBS - see note on page 6**



Tooth measurements are within plus or minus 5% of stated value.



# ORDER FORM

## REGAL-D-BLEND® PRODUCTS:



Order No.: \_\_\_\_\_ Date: \_\_\_\_\_  
 Bill To: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_  
 Phone # \_\_\_\_\_ Fax # \_\_\_\_\_  
 Ship To: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

**Item No.**  
 079-50-0901 Regal-D-Blend Anterior Teeth (1 x 6)  
 084-50-0903 Regal-D-Blend Maxillary Centrals Shade Guide in 6 Shades  
 084-50-0904 Regal-D-Blend Anterior Working Mold Guide  
 Contains one set of each mold - 35 Anteriors, all one shade  
 079-58-0003 Justi Blend Posterior Teeth (1 x 8)  
 084-59-0011 Justi Blend Posterior Working Mold Guide  
 Contains one set of each mold - 36 posteriors, all one shade

1 x 6 27 UPPER ANTERIOR REGAL-D-BLEND® MOLDS  Boxes

Shades	AX/200	AA/202	BB/204	CC/206	DD/208	EE/210
<b>Molds</b>						
A4						
A7						
B3						
B8						
C1						
C2						
C5						
C6						
D1						
D2						
D3						
D6						
E1						
E2						
E3						
E6						
E8						
F1						
F4						
G5						
G7						
H2						
H4						
J2						
J3						
J5						
K4						
Total						

1 x 6 8 LOWER ANTERIOR REGAL-D-BLEND® MOLDS  Boxes

Shades	AX/200	AA/202	BB/204	CC/206	DD/208	EE/210
<b>Molds</b>						
1						
2						
3						
4						
5						
6						
7						
8						
Total						

Anterior Cards (1 x 6) - 14 cards per box.

Total Upper Anteriors 1 x 6	
Total Lower Anteriors 1 x 6	
Total Anteriors 1 x 6	

Posterior Cards (1 x 8) - 20 cards per box.

Total Upper Posteriors 1 x 8	
Total Lower Posteriors 1 x 8	
Total Posteriors 1 x 8	

Item #079-50-0901

Item #079-58-0003

Remarks \_\_\_\_\_

Signed \_\_\_\_\_

1 x 8 18 UPPER POSTERIOR JUSTI BLEND® MOLDS  Boxes

Shades	AX/200	AA/202	BB/204	CC/206	DD/208	EE/210
<b>Molds</b>						
0°						
29MF						
31MF						
33MF						
10°						
30F						
32F						
34F						
20°						
29L						
31L						
33L						
29M						
31M						
33M						
ANATOMICAL MEDIUM - 33°						
30M						
32M						
34M						
ANATOMICAL LONG - 33°						
30L						
32L						
34L						
Total						

1 x 8 18 LOWER POSTERIOR JUSTI BLEND® MOLDS  Boxes

Shades	AX/200	AA/202	BB/204	CC/206	DD/208	EE/210
<b>Molds</b>						
0°						
29MF						
31MF						
33MF						
10°						
30F						
32F						
34F						
20°						
29L						
31L						
33L						
29M						
31M						
33M						
ANATOMICAL MEDIUM - 33°						
30M						
32M						
34M						
ANATOMICAL LONG - 33°						
30L						
32L						
34L						
Total						

NOTE: Posteriors are also available in shades 212, 214 and 216 by special request.

This is a sample of our tooth order form. You may use a photo copy of this form for your first order.

See page 6 for details on Justi Blend posterior teeth.

# Regal-D-Blend® Articulation Chart

Anteriors		Posteriors			
Upper Molds	Articulates Lower Anteriors	Zero Degree MF	10 Degree F	20 Degree	Anatomical 33 Degree Short, Med., Long
A4	4	31	32	29L	30
A7	7	31	-	33M	32
B3	3	31	32	31M	30
B8	8	31	-	33M	32
C1	1	29	30	29L	30
C2	2	29	30	29L	30
C5	5	31	32	31L	32
C6	6	31	32	31L	32
D1	1	29	30	29M	30
D2	2	31	32	31L	30
D3	3	31	32	29L	30
D6	6	33	34	33L	34
E1	1	29	30	29L	30
E2	2	31	30	29L	30
E3	3	31	32	31L	30
E6	6	31	32	31L	32
E8	8	33	34	33L	32
F1	1	29	30	29L	30
F4	4	31	-	31M	30
G5	5	31	32	31L	32
G7	7	31	32	31L	32
H2	2	31	30	29L	30
H4	4	31	32	31L	30
J2	2	29	32	29L	30
J3	3	29	32	31L	30
J5	5	31	34	33L	32
K4	4	31	-	29L	30



Item# 084-50-0904  
Regal-D-Blend Anterior Working Mold Guide.  
One set of each mold - 35 anteriors, all one shade.

Item# 084-50-0903  
Regal-D-Blend Maxillary Centrals Shade Guide in 6 shades.

## \* Posterior Teeth

Regal-D-Blend® Anterior Hardened Plastic Teeth have a perfect mate for the posterior arrangement with Justi Blend® Posterior Teeth. The posteriors pictured in this mold chart are the Justi Blend® Posteriors. Justi Blend® Posteriors, like Regal-D-Blend® Anteriors, meet and exceed all ADA Accepted Specifications.

## Anterior Mold Comparison With Dentsply's Bioblend®, Trublend®, Portrait®

Regal-D-Blend®	Bioblend® IPN® Bioblend® Trublend®, Portrait®	Regal-D-Blend®	Bioblend® IPN® Bioblend® Trublend®, Portrait®
<b>Uppers</b>		<b>Uppers</b>	
A4	12E	F4	43F
A7	12G	G5	45F
B3	32E	G7	45H
B8	32G	H2	55D
C1	21C	H4	55F
C2	21D	J2	62D
C5	21E	J3	55D
C6	21X	J5	62G
D1	22C	K4	75E
D2	13E	<b>Lowers</b>	
D3	22E	1	C
D6	22H	2	F
E1	42C	3	H
E2	42D	4	N
E3	42F	5	O
E6	42G	6	P
E8	42H	7	R
F1	43D	8	S

## Shade Comparison with Major Brands

Color	Regal-D-Blend®	Justi Blend®	Bio-blend®	Bio-form®	Tru-Blend™	Vita®	Por-trait™	Ivoclar Vivadent
White/Gray	AX	200	100	B59	T2	A-1	P-1	-
Yellow/Pink	AA	202	102	B53	T3	A2-B2	P-2	1A
Yellow/Pink	BB	204	104	B92	T4	A2	P-2	1A, 2A
Gray/Pink	CC	206	106	B69	T5	C2-D3	P69	-
Yellow/Brown	DD	208	108	B62	T12	A3-D3	P-3	1C
Gray	EE	210	110	B66/93	T15	A35	P33/P66	1C

This shade comparison chart is to be used in selecting approximate color matches. Slight variations occur due to the difference in the texture and character of the respective materials. This chart should be used as an aid to the skill and preferences of the dentist or technician who must make the final shade choice.

**Regal-D-Blend® Anteriors match the central, lateral and cuspid, right & left respectively. Each central, lateral and cuspid are the same shade in three different tones thereby creating a natural blended effect.**

**Justi Blend® Posteriors match the right or left cuspid body of Regal-D-Blend in tone and shade.**

**Regal-D-Blend® teeth meet and exceed ADA Specification 15 for Certification.**

Distributed by:

## Instructions for use of synthetic resin teeth

Note: If you would like to receive written instructions for obtaining optimum results with our teeth please contact:  
Justi Products, A Division of American Tooth Industries,  
1200 Stellar Drive, Oxnard, CA 93033-2404  
Toll Free (800) 235-4639, California (805) 487-9868  
Fax (805) 483-8482 • Website: [www.americantooth.com](http://www.americantooth.com)

MANUFACTURED BY:



**Justi Products**

A Brand Division of American Tooth Industries

1200 Stellar Drive, Oxnard, CA 93033-2404

(805) 487-9868 • FAX (805) 483-8482

E-mail: [info@americantooth.com](mailto:info@americantooth.com) • Website: [www.americantooth.com](http://www.americantooth.com)