



DEDICATED TO OUR TOOTH EXPRESS® CUSTOMERS AND TO OUR INDEPENDENT TOOTH EXPRESS® DEALER/LICENSEES

TOOTH EXPRESS®

Published by American Tooth Industries to promote Justi, Pozzi, LeoneAmerica, Major America and Regal Dental Products to Dental Specialists throughout the world.
 Tooth Express® • 1200 Stellar Drive • Oxnard, California 93033-2404 •
 Toll Free # (California and National) (800) 628-1437 • FAX (805) 483-8482
 Web site: www.americantooth.com • E-mail: info@americantooth.com

News




THE **NAKED** TRUTH...

JUSTI EDUCATION: FOR TECHNICIANS & DENTISTS - TODAY & TOMORROW

At American Tooth Industries, we know that the dental industry is very demanding in terms of the knowledge required to enter the field, and the continuing education required to keep up with the ever-changing technology.

Our Justi line of educational materials is designed to get you there, and to keep you there:



1 **P.R.E.C.D.T.™**

A comprehensive digitized tutorial to help technicians prepare for the CDT Exams.

2 **TECHNICIAN'S MANUALS**

Comprehensive manuals, prepared by U.S. Air Force Dental Laboratory Instructors, were written for technicians who are studying various fabrication procedures used in the dental lab as well as the basic sciences needed for successful communication with dentists. They are also recommended for technicians who are preparing to take certification exams.

A must for professional dental technicians!



3 **PAPILLA GAUGE - Two Dimensional Gauge**

The position of the papilla can provide great diagnostic value in reproducing tooth position in the construction of a denture. The Justi Papilla Gauge will accurately measure the horizontal and vertical plane to define the position of the incisal edges and the central incisors.



4 **ARTIFICIAL MOLD SELECTOR**

A unique tooth selecting system, designed to aid the dental and laboratory professional in the selection of specific Justi tooth molds. It provides 100% assurance in accuracy of mold selection. Designed and patented by Justi.



5 **JUSTI POCKET LENS**



See all the details around you without having to worry about misplacing an expensive magnifying glass.

6 **TECHNICIAN'S UTILITY RULER**

An invaluable aid to the lab technician. Fabrication of occlusal rims, verification of the maxillary central incisor position, determination of normal or cross-bite situations, determination of posterior molds, measuring for the six anterior teeth and related information.



7 **UTILITY TEMPLATE**

The Justi Utility Template is a quality stainless occlusal plane guide tool specifically designed for the needs of today's productive technician. Perfect for tooth set-up and occlusal bite rim fabrication. A guide relation of 80 degree angled sides are calibrated with a millimeter rule guide measuring 0 to 70 mm with markings at 18 and 22 mm.

8 **JUSTI CROWN GAUGE**

A traditional stainless crown caliper which can also be used for measuring the buccal flange thickness on an acrylic denture base. Measurements in 1/10 mm increments up to 10 mm.



9 **JUSTI CALIPER**

A compass-style measuring device used for measuring parameter landmarks on a working model or facial dimensions to aid in determining tooth mold selection.



MORE EDUCATION...



NEW FROM...



ANATOMICAL TOOTH MODEL 10

A large resin, two-rooted molar, approximately 8" tall. Depicts decay into enamel, dentin, and pulp. The beginning and progression of decay resulting in an abscess. Includes a removable crown which reveals a molar crown prep.



ANATOMICAL STUDY MODELS 11

The perfect tool for patient education. The crystal clear model shows each tooth in its entirety as well as its relationship to the rest of the dentition. Each individual tooth may be removed to explain treatment. The patented Vertex hinge guarantees a lifetime of use with these models.



The anatomical Orthodontic Study Model was made showing imperfect dentition, allowing the practitioner to show the patient some of the common recommendations for treatment. *Soft tissues on both models.*

DENTERMS - DENTAL DICTIONARY 12

A glossary of dental terminology perfect for the prosthodontic dental team. Much of the "language of dentistry" is defined. For use by schools of dental laboratory technology, schools of dental assisting, dental hygiene and anyone who is interested in dental health and is seeking a basic knowledge of dental anatomy.



JUSTI TOOTH MANUAL 13

A manual for the selection, articulation and fabrication of dentures using Justi plastic teeth. A valuable guide for the dental laboratory technician. State-of-the-art illustrations and photographs describing some of the types of removable prosthodontic appliances and techniques used in dentistry today. A valuable aid to both teaching and learning.

...AND MANY MORE



"Believe Me!
The Best
for Labs..."
Bruno Pozzi

The new Pozzi Malibu series of molar bands are fully proportional in all dimensions: diameter, height and thickness. This provides maximum strength and prevents cement washout due to gapping.

Pozzi Malibu bands are designed with a regular temper to allow for easier "stretch" adaptation. They are available in separate lefts and rights for better seating and a closer adaptation to the natural anatomy.

Pozzi Malibu bands utilize Unitek® standard sizing, making switching easy and painless. If you are tired of paying high prices for your bands, then the Pozzi Malibu line is perfect for you!

EXPANSOR ASSORTMENT KIT FOR REMOVABLE ORTHODONTIC APPLIANCES

Contains 270 assorted pieces from the 500, 800 and 900 series of removable Expandors. It is perfect as a high-powered starter kit, economically priced, for the construction of any type of removable appliance which requires an Expansor that is embedded in acrylic.

The heavy duty storage tray is perfect for maintaining an accurate inventory of Expandors.

Full Assortment at Factory Prices






NEW EXPANSOR™ ATLAS FROM POZZI & LEONE

Call to Request Your Copy 800 628-1437

The new Expansor Series Atlas for 2002 features a complete selection of Expansors for both fixed and removable orthodontic appliances, as well as all of the accessories require for fabrication.

EXPANSORS FOR REMOVABLE APPLIANCES



Three Directional Expansor - 14mm



Fan-Type Rapid Expansor - 9mm

ACCESSORIES

EXPANSORS FOR FIXED APPLIANCES



Pozzi Ortho Colors




Pozzi Muscle U

ORMALAB!

ORMALAB 95

When extreme hardness is required. For use with self curing and heat curing acrylic resins.



Call for Special Introductory Factory Price

These polysiloxane based, condensation curing materials were designed exclusively for extra-oral use in the dental laboratory.

Ormalab 95 is indispensable whenever hardness and heat-resistance are required. It is ideal for matrices for fixed and total prosthesis, duplications, prosthesis repair, and matrices for temporary crowns and bridges. It is extremely heat resistant - over 120° C without distortion. It is extremely hard, comparable to gypsum during the milling phase.

Ormalab 75 is intended for matrices for fixed and total prosthesis, prosthesis repair, and matrices for temporary crowns and bridges. It is compatible with cold curing acrylic resins, and features long lasting dimensional stability and high compressive strength and recovery from deformation.

ORMALAB 75

When flexibility is the primary concern. For use with cold curing acrylic resins.





Great New Products & Important Information!

For the Dealer or
Tooth Express® Store
nearest you, please call us at
(805) 487-9868 or
Toll Free (800) 628-1437.



PRSR STD
MAIL
U.S. POSTAGE
PAID
OXNARD, CA
PERMIT NO. 539

TOOTH EXPRESS®
Published by:
AMERICAN TOOTH INDUSTRIES
1200 Stellar Drive
Oxnard, CA USA
93033-2404